


Anti-Phospho-C-Jun (T91) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P05412

Item Number: CM0001161



Introduction

Rabbit monoclonal antibody against phospho-c-Jun (T91) for WB, IHC-P, ELISA in human and mouse. Targets transcription factor Jun, a component of AP-1 complex, with roles in cell proliferation, apoptosis, and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Jun-T91 Rabbit mAb
Remarks/Alias	AP1; p39; AP-1; cJUN; c-Jun; Phospho-c-Jun-T91
Species	Human
GeneID (Human)	3725
GeneID	3725
Immunogen	A synthetic phosphorylated peptide around T91 of human c-Jun (P05412)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05412
Protein Weight	36kDa
Shipping	Ice bag