


Anti-Phospho-C-Myc-S373 Rabbit Polyclonal Antibody For Wb - P01106

Item Number: CM0001565



Introduction

Rabbit polyclonal antibody specific for human c-Myc phosphorylated at Ser373. Suitable for Western blot analysis. Targets the proto-oncogene transcription factor c-Myc (~51 kDa), involved in cell growth and oncogenesis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Myc-S373 Rabbit pAb
Remarks/Alias	MRTL; MYCC; c-Myc; bHLHe39; Phospho-c-Myc-S373
Species	Human
Gene ID (Human)	4609
Gene ID	4609
Immunogen	A phospho specific peptide corresponding to residues surrounding S373 of human c-Myc
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P01106
Protein Weight	51kDa
Shipping	Ice bag