


Anti-Phospho-C-Myc (S62) Rabbit Polyclonal Antibody For Wb, Ihc-P, And Elisa - P01106

Item Number: CM0001103



Introduction

Phospho-specific rabbit polyclonal antibody targeting c-Myc phosphorylated at Ser62. Validated for WB, IHC-P, and ELISA in human, mouse, and rat. Ideal for cancer research and signal transduction studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Myc-S62 Rabbit pAb
Remarks/Alias	MRTL; MYCC; c-Myc; bHLHe39; Phospho-c-Myc-S62
Species	Human
GeneID (Human)	4609
Gene ID	4609
Immunogen	Synthetic peptide A phospho specific peptide corresponding to residues surrounding S62 of human MYC (NP_001341799.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01106
Protein Weight	51kDa
Shipping	Ice bag