


Anti-Phospho-C-Myc (T358) Polyclonal Antibody For Wb, Elisa - P01106

Item Number: CM0001406



Introduction

Rabbit polyclonal antibody against phosphorylated c-Myc at Thr358. Suitable for Western blot and ELISA to detect human c-Myc phosphorylation, a key transcription factor and oncogene involved in cell growth and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-c-Myc-T358 Rabbit pAb
Remarks/Alias	MRTL; MYCC; c-Myc; bHLHe39; Phospho-c-Myc-T358
Species	Human
GeneID (Human)	4609
GeneID	4609
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding T358 of human c-Myc
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01106
Protein Weight	51kDa
Shipping	Ice bag