


Anti-Phospho-C-Myc (T58) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P01106

Item Number: CM0001544



Introduction

Rabbit polyclonal antibody targeting phosphorylated c-Myc (T58), validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for Myc signaling, cell cycle, and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-c-Myc-T58 Rabbit pAb
Remarks/Alias	MRTL; MYCC; c-Myc; bHLHe39; Phospho-c-Myc-T58
Species	Human
Gene ID (Human)	4609
Gene ID	4609
Immunogen	A synthetic phosphorylated peptide around T58 of mouse MYC (NP_001170823.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01106
Protein Weight	51kDa
Shipping	Ice bag