


Anti-Phospho-Camk2-T286 Polyclonal Antibody For Wb, Elisa - Q13554 / Q13557 / Q13555 / Q9Uqm7

Item Number: CM0001858



Introduction

Rabbit polyclonal antibody specifically recognizing phosphorylated T286 of CaMKII alpha (Camk2a) with cross-reactivity to human, mouse, and rat. Validated for WB and ELISA, ideal for studying synaptic plasticity and CaMKII activation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Camk2-T286 Rabbit pAb
Remarks/Alias	CAMKA; MRD53; MRT63; CaMKIIalpha; CaMKIINalpha; CAM2; CAMK2; CAMKB; MRD54; CaMKIibeta; CAMKD; CAMK; CAMKG; MRD59; CAMK-II; Phospho-Camk2-T286
Species	Human
Gene ID (Human)	815
Gene ID	816/817/818/815
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T286 of human Camk2a (NP_741960.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13554/Q13557/Q13555/Q9UQM7
Protein Weight	54kDa
Shipping	Ice bag