


Anti-Phospho-Caveolin-1-Y14 Rabbit Polyclonal Antibody For Wb And Elisa - Q03135

Item Number: CM0001984



Introduction

Rabbit polyclonal antibody specific for phosphorylated Caveolin-1 at tyrosine 14. Validated for Western Blot and ELISA, with cross-reactivity to human and mouse. Ideal for caveolae signaling and T-cell activation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Caveolin-1-Y14 Rabbit pAb
Remarks/Alias	CGL3; PPH3; BSCL3; LCCNS; VIP21; MSTP085; Phospho-Caveolin-1-Y14
Species	Human
GeneID (Human)	857
GeneID	857
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around Y14 of human Caveolin-1 (NP_001744.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q03135
Protein Weight	20kDa
Shipping	Ice bag