


Anti-Phospho-Cbl-Y774 Rabbit Pab - P22681

Item Number: CM0000832



Introduction

Phospho-CBL-Y774 Rabbit Polyclonal Antibody validated for WB and ELISA, specifically detects phosphorylated CBL at tyrosine 774 in human samples. Ideal for cancer and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CBL-Y774 Rabbit pAb
Remarks/Alias	CBL2; NSLL; C-CBL; RNF55; FRA11B; Phospho-CBL-Y774
Species	Human
GeneID (human)	867
GeneID	867
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y774 of human CBL (NP_005179.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P22681
Protein Weight	100kDa
Shipping	Ice bag