


Anti-Phospho-Cd115/Csf-1R-Y723 Rabbit Polyclonal Antibody For Wb And Elisa - P07333

Item Number: CM0001443



Introduction

Phospho-specific rabbit polyclonal antibody targeting human CSF1R (CD115) phosphorylated at Y723. Validated for WB and ELISA. Ideal for investigating CSF1/IL34 signaling, macrophage biology, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CD115/CSF-1R-Y723 Rabbit pAb
Remarks/Alias	FMS; CSFR; FIM2; GPSC; HDLS; C-FMS; CD115; HDLS1; CSF-1R; BANDDOS; M-CSF-R; Phospho-CD115/CSF-1R-Y723
Species	Human
GeneID (Human)	1436
GeneID	1436
Immunogen	Synthetic peptide (a synthetic phosphorylated peptide around Y723 of human CSF1R, NP_005202.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P07333
Protein Weight	108kDa
Shipping	Ice bag