

Anti-Phospho-Cd117/C-Kit-Y721 Rabbit Pab For Wb - P10721

Item Number: CM0000449



Introduction

Anti-Phospho-CD117/c-Kit-Y721 Rabbit Polyclonal Antibody, validated for Western Blot in Human. Detects phosphorylated tyrosine 721 of KIT, a receptor tyrosine kinase crucial for cell survival, proliferation, and hematopoiesis. Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CD117/c-Kit-Y721 Rabbit pAb
Remarks/Alias	PBT; SCFR; C-Kit; CD117; MASTC; Phospho-CD117/c-Kit-Y721
Species	Human
GeneID (Human)	3815
GeneID	3815
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding Y721 of human KIT
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P10721
Protein Weight	110kDa
Shipping	Ice bag