

Anti-Phospho-Cdc25C (S216) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P30307

Item Number: CM0001157



Introduction

Phospho-specific cdc25C (S216) rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat targets. Ideal for studying cell cycle progression and mitotic regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-cdc25C-S216 Rabbit pAb
Remarks/Alias	CDC25; PPP1R60; Phospho-cdc25C-S216
Species	Human
GeneID (Human)	995
GeneID	995
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S216 of human cdc25C (NP_001781.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30307
Protein Weight	53kDa
Shipping	Ice bag