


Anti-Phospho-Cdc37-S13 Rabbit Monoclonal Antibody For Wb, Elisa - Q16543

Item Number: CM0001589



Introduction

Phospho-specific rabbit monoclonal antibody targeting CDC37 phosphorylated at Ser13. Validated for Western blot and ELISA in human, mouse, and rat samples. Ideal for studying Hsp90 co-chaperone function and kinase stabilization pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDC37-S13 Rabbit mAb
Remarks/Alias	P50CDC37; Phospho-CDC37-S13
Species	Human
GeneID (Human)	11140
GeneID	11140
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16543
Protein Weight	44kDa
Shipping	Ice bag