


Anti-Phospho-Cdk1-T14 Rabbit Polyclonal Antibody For Wb, Elisa - P06493

Item Number: CM0000610



Introduction

Rabbit polyclonal antibody recognizing phospho-Thr14 of human CDK1 (CDC2), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for cell cycle, mitosis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK1-T14 Rabbit pAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; Phospho-CDK1-T14
Species	Human
GeneID (Human)	983
GeneID	983
Immunogen	A synthetic phosphorylated peptide around T14 of human CDK1 (NP_001777.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06493
Protein Weight	34kDa
Shipping	Ice bag