


Anti-Phospho-Cdk1-T161+Cdk2/Cdk3-T160 Rabbit Monoclonal Antibody For Wb - P06493

Item Number: CM0001838



Introduction

Rabbit monoclonal antibody recognizing phosphorylated CDK1-T161, CDK2-T160, and CDK3-T160. Suitable for Western blot, IF-P, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for cell cycle, mitosis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK1-T161+CDK2/CDK3-T160 Rabbit mAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; Phospho-CDK1-T14
Species	Human
GeneID (Human)	983
GeneID	983
Immunogen	A synthetic phosphorylated peptide around T161&T160 of human CDK1 (NP_001777.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06493
Protein Weight	34 kDa
Shipping	Ice bag