


Anti-Phospho-Cdk1-Y15 Polyclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - P06493

Item Number: CM0001850



Introduction

Phospho-CDK1-Y15 Rabbit Polyclonal Antibody validated for WB, IP, IF-P, IHC-P, ELISA in human, mouse, rat. Recognizes CDK1 phosphorylated at Tyr15, a key cell cycle regulator. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CDK1-Y15 Rabbit pAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; Phospho-CDK1-Y15
Species	Human
GeneID (Human)	983
GeneID	983
Immunogen	Synthetic phosphorylated peptide around Y15 of human CDK1 (NP_001777.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06493
Protein Weight	34kDa
Shipping	Ice bag