


## Anti-Phospho-Cdk9 (T186) Rabbit Polyclonal Antibody For Wb, Elisa - P50750

Item Number: CM0000615



### Introduction

Phospho-CDK9 (T186) rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat CDK9. Ideal for transcription regulation and kinase signaling research.

[Learn More](#)

| Parameter        | Value                                                                                          |
|------------------|------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-CDK9-T186 Rabbit pAb                                                                   |
| Remarks/Alias    | TAK; C-2k; CTK1; CDC2L4; PITALRE; Phospho-CDK9-T186                                            |
| Species          | Human                                                                                          |
| GeneID (Human)   | 1025                                                                                           |
| GeneID           | 1025                                                                                           |
| Immunogen        | Synthetic Peptide: A synthetic phosphorylated peptide around T186 of human CDK9 (NP_001252.1). |
| Source           | Rabbit                                                                                         |
| Category         | Phosphorylated Antibodies                                                                      |
| Application      | WB, ELISA                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                              |
| SWISS            | P50750                                                                                         |
| Protein Weight   | 43kDa                                                                                          |
| Shipping         | Ice bag                                                                                        |