

Anti-Phospho-Cfl1-Y139 Rabbit Polyclonal Antibody For Ihc-P, Elisa - P23528

Item Number: CM0001533



Introduction

Phospho-CFL1-Y139 Rabbit pAb from CamelBio detects phosphorylated cofilin-1 (Y139) in human, mouse, and rat for IHC-P and ELISA. Ideal for actin dynamics and cell signaling research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-CFL1-Y139 Rabbit pAb
Remarks/Alias	CFL; cofilin; HEL-S-15; Phospho-CFL1-Y139
Species	Human
GeneID (Human)	1072
GeneID	1072
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23528
Protein Weight	19kDa
Shipping	Ice bag