


Anti-Phospho-Chk1 (S280) Monoclonal Antibody For Wb, Elisa - O14757

Item Number: CM0001670



Introduction

Phospho-specific rabbit monoclonal antibody detecting human Chk1 phosphorylated at serine 280. Suitable for Western blot and ELISA applications. Designed for cell cycle checkpoint and DNA damage response research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Chk1-S280 Rabbit mAb
Remarks/Alias	CHK1; Phospho-Chk1-S280
Species	Human
Gene ID (Human)	1111
Immunogen	A phospho synthetic peptide corresponding to residues surrounding S280 of Human Chk1.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag