


Anti-Phospho-Chk1-S280 Rabbit Monoclonal Antibody For Wb, Elisa - O14757

Item Number: CM0001468



Introduction

Phospho-Chk1-S280 rabbit monoclonal antibody validated for Western blot (WB) and ELISA. Detects human Chk1 phosphorylated at serine 280. Suitable for studying DNA damage response, cell cycle checkpoints, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Chk1-S280 Rabbit mAb
Remarks/Alias	CHK1; Phospho-Chk1-S280
Species	Human
GeneID (Human)	1111
GeneID	1111
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag