


Anti-Phospho-Chk1-S296 Rabbit Monoclonal Antibody For Wb, Elisa - O14757

Item Number: CM0000717



Introduction

High-quality rabbit monoclonal antibody against phospho-Chk1-S296 (human) for Western Blot and ELISA. Essential for studying DNA damage response and checkpoint signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Chk1-S296 Rabbit mAb
Remarks/Alias	CHK1; Phospho-Chk1-S296
Species	Human
GeneID (human)	1111
GeneID	1111
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S296 of human Chk1 (O14757).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag