


## Anti-Phospho-Chk1 (S317) Rabbit Monoclonal Antibody For Wb, Elisa - O14757

Item Number: CM0001461



### Introduction

Rabbit monoclonal antibody targeting phospho-Chk1 (Ser317). Validated for WB and ELISA, it enables detection of activated Chk1 in human samples for DNA damage and cell cycle checkpoint research.

[Learn More](#)

| Parameter        | Specification                                                              |
|------------------|----------------------------------------------------------------------------|
| Product Name     | Phospho-Chk1-S317 Rabbit mAb                                               |
| Remarks/Alias    | CHK1; Phospho-Chk1-S317                                                    |
| Species          | Human                                                                      |
| Gene ID (Human)  | 1111                                                                       |
| Gene ID          | 1111                                                                       |
| Immunogen        | A synthetic phosphorylated peptide around S317 of human Chk1 (NP_001265.2) |
| Source           | Rabbit                                                                     |
| Category         | Phosphorylated Antibodies                                                  |
| Application      | WB, ELISA                                                                  |
| Cross Reactivity | Human                                                                      |
| SWISS            | O14757                                                                     |
| Protein Weight   | 54kDa                                                                      |
| Shipping         | Ice bag                                                                    |