


Anti-Phospho-Chk1-S345 Rabbit Polyclonal Antibody For Wb And Elisa - O14757

Item Number: CM0000661



Introduction

Phospho-specific rabbit polyclonal antibody targeting human Chk1 at Ser-345. Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Involved in DNA damage checkpoint and cell cycle arrest pathways.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-Chk1-S345 Rabbit pAb
Remarks/Alias	CHK1; Phospho-Chk1-S345
Species	Human
Gene ID (Human)	1111
Gene ID	1111
Immunogen	A synthetic phosphorylated peptide around S345 of human Chk1 (NP_001265.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag