

Anti-Phospho-Chk2-T68 Rabbit Polyclonal Antibody For Wb, Elisa - O96017

Item Number: CM0000772



Introduction

Rabbit polyclonal antibody specific for phospho-T68 of human Chk2 (CHEK2), validated for WB and ELISA. Ideal for DNA damage signaling and cell cycle checkpoint research.

[Learn More](#)

Parameter	Details
Product Name	Phospho-Chk2-T68 Rabbit pAb
Remarks/Alias	CDS1; CHK2; LFS2; RAD53; hCds1; HuCds1; PP1425; Phospho-Chk2-T68
Species	Human
Gene ID	11200
Immunogen	Synthetic phosphorylated peptide around T68 of human Chk2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O96017
Protein Weight	61kDa
Shipping	Ice bag