

Anti-Phospho-Creb1-S133 Rabbit Monoclonal Antibody For IHC-P, Elisa - P16220

Item Number: CM0001240



Introduction

This phospho-CREB1 (Ser133) antibody detects human, mouse, and rat proteins. Ideal for IHC-P and ELISA, it is a valuable tool for studying transcription regulation, circadian rhythms, and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CREB1-S133 Rabbit mAb
Remarks/Alias	CREB; CREB-1; Phospho-CREB1-S133
Species	Human
GeneID	1385
Immunogen	Synthetic peptide (A synthetic phosphorylated peptide around Ser133 of human CREB1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16220
Protein Weight	35kDa
Shipping	Ice bag