

Anti-Phospho-Creb1-S133 Rabbit Pab For Wb, Ihc-P, Elisa - P16220

Item Number: CM0001169



Introduction

Rabbit polyclonal antibody against phosphorylated CREB1 (Ser133). Validated for WB, IHC-P, ELISA with reactivity across human, mouse, and rat. Ideal for transcription, circadian rhythm, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CREB1-S133 Rabbit pAb
Remarks/Alias	CREB; CREB-1; Phospho-CREB1-S133
Species	Human
GeneID (human)	1385
GeneID	1385
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S133 of human CREB1 (NP_604391.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16220
Protein Weight	35kDa
Shipping	Ice bag