

Anti-Phospho-Creb1 (Ser129) Rabbit Polyclonal Antibody For Wb, Elisa - P16220

Item Number: CM0000645



Introduction

Phospho-specific rabbit polyclonal antibody targeting CREB1 phosphorylated at Ser129. Validated for WB and ELISA in human, mouse, and rat samples. Suitable for studying cAMP-mediated transcription and circadian rhythm regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CREB1-S129 Rabbit pAb
Remarks/Alias	CREB; CREB-1; Phospho-CREB1-S129
Species	Human
Gene ID (Human)	1385
Gene ID	1385
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S129 of human CREB1 (NP_604391.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16220
Protein Weight	35kDa
Shipping	Ice bag