


Anti-Phospho-Crk/Crkii (Y221) Rabbit Polyclonal Antibody For Wb And Elisa - P46108

Item Number: CM0001788



Introduction

Rabbit polyclonal antibody specific to phosphorylated Y221 of human CRK/CRKII. Validated for WB and ELISA; reacts with human, mouse, rat. Ideal for studying CRK-mediated cell adhesion and migration signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CRK/CRKII-Y221 Rabbit pAb
Remarks/Alias	p38; CRKII; Phospho-CRK/CRKII-Y221
Species	Human
GeneID	1398
Immunogen	A synthetic phosphorylated peptide around Y221 of human CRK/CRKII (NP_058431.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46108
Protein Weight	34kDa
Shipping	Ice bag