

Anti-Phospho-Crkl (Y207) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P46109

Item Number: CM0001679



Introduction

Phospho-CRKL-Y207 Rabbit pAb validated for WB, IHC-P, ELISA in human, mouse, rat. A synthetic phosphopeptide around Y207 of CRKL (P46109) was used as immunogen. Suitable for studying intracellular signal transduction.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CRKL-Y207 Rabbit pAb
Remarks/Alias	CRKL; Phospho-CRKL-Y207
Species	Human
GeneID (Human)	1399
GeneID	1399
Immunogen	A synthetic phosphorylated peptide around Y207 of human CRKL (NP_005198.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46109
Protein Weight	34kDa
Shipping	Ice bag