


Anti-Phospho-Crmp2/Dpysl2-T514 Rabbit Polyclonal Antibody For Wb, Elisa - Q16555

Item Number: CM0000721



Introduction

Phospho-CRMP2/DPYSL2-T514 Rabbit pAb validated for WB, ELISA; targets phosphorylated T514 of human CRMP2. Cross-reacts with mouse and rat. Ideal for neuronal development, axon guidance, and semaphorin signaling studies.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-CRMP2/DPYSL2-T514 Rabbit pAb
Remarks/Alias	DRP2; N2A3; CRMP2; DRP-2; ULIP2; CRMP-2; DHPRP2; ULIP-2; Phospho-CRMP2/DPYSL2-T514
Species	Human
GeneID (Human)	1808
GeneID	1808
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around T514 of human CRMP2/DPYSL2 (NP_001377.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16555
Protein Weight	62kDa
Shipping	Ice bag