


Anti-Phospho-Csf-1R (Y708) Rabbit Polyclonal Antibody For Wb, Elisa - P07333

Item Number: CM0001578



Introduction

Rabbit polyclonal antibody specific for phosphorylated CSF-1R at Y708. Validated for Western blot and ELISA on human samples. Ideal for studying macrophage signaling, innate immunity, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-CD115/CSF-1R-Y708 Rabbit pAb
Remarks/Alias	FMS; CSFR; FIM2; GPSC; HDLS; C-FMS; CD115; HDLS1; CSF-1R; BANDDOS; M-CSF-R; Phospho-CD115/CSF-1R-Y708
Species	Human
GeneID (Human)	1436
GeneID	1436
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y708 of human CSF1R (NP_005202.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P07333
Protein Weight	108kDa
Shipping	Ice bag