


## Anti-Phospho-Cytokeratin 8 (Krt8)-S74 Rabbit Polyclonal Antibody For Wb And Elisa - P05787

Item Number: CM0001639



### Introduction

Phospho-specific rabbit polyclonal antibody detects human Keratin 8 phosphorylated at Serine 74. Validated for WB and ELISA, suitable for studying intermediate filament dynamics and epithelial integrity.

[Learn More](#)

| Parameter        | Value                                                                               |
|------------------|-------------------------------------------------------------------------------------|
| Product Name     | Phospho-Cytokeratin 8 (KRT8)-S74 Rabbit pAb                                         |
| Remarks/Alias    | K8; KO; CK8; CK-8; CYK8; K2C8; CARD2; Phospho-Cytokeratin 8 (KRT8)-S74              |
| Species          | Human                                                                               |
| GeneID (Human)   | 3856                                                                                |
| GeneID           | 3856                                                                                |
| Immunogen        | A synthetic phosphorylated peptide around S74 of human Cytokeratin 8 (NP_002264.1). |
| Source           | Rabbit                                                                              |
| Category         | Phosphorylated Antibodies                                                           |
| Application      | WB, ELISA                                                                           |
| Cross Reactivity | Human                                                                               |
| SWISS            | P05787                                                                              |
| Protein Weight   | 54 kDa/56 kDa                                                                       |
| Shipping         | Ice bag                                                                             |