

Anti-Phospho-Ddr1 (Y792) Rabbit Polyclonal Antibody For Wb, Elisa - Q08345

Item Number: CM0000869



Introduction

Phospho-DDR1 (Y792) rabbit polyclonal antibody validated for WB and ELISA. Detects human, mouse, rat DDR1 tyrosine kinase (101 kDa). Essential for collagen signaling and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-DDR1-Y792 Rabbit pAb
Remarks/Alias	CAK, DDR, NEP, HGK2, PTK3, RTK6, TRKE, CD167, EDDR1, MCK10, NTRK4, PTK3A, Phospho-DDR1-Y792
Species	Human
Gene ID (Human)	780
Gene ID	780
Immunogen	Synthetic phosphorylated peptide around Y792 of human DDR1 (NP_054699.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08345
Protein Weight	101 kDa
Shipping	Ice bag