

Anti-Phospho-Dnmt1 (S714) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P26358

Item Number: CM0000479



Introduction

Phospho-specific DNMT1 (Ser714) rabbit polyclonal antibody validated for Western blot, IF/ICC, and ELISA. Detects human DNMT1 phosphorylated at Ser714. Suitable for epigenetic inheritance, gene silencing, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-DNMT1-S714 Rabbit pAb
Remarks/Alias	AIM; DNMT; MCMT; CXXC9; HSN1E; ADCADN; m.Hsal; Phospho-DNMT1-S714
Species	Human
Gene ID (Human)	1786
Gene ID	1786
Immunogen	A synthetic phosphorylated peptide around S714 of human DNMT1 (NP_001370.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P26358
Protein Weight	183kDa
Shipping	Ice bag