


Anti-Phospho-Drp1-S616 Polyclonal Antibody For Wb And Elisa - O00429

Item Number: CM0000719



Introduction

Phospho-DRP1-S616 Rabbit pAb specifically detects human DNM1L phosphorylated at Ser616. Validated for WB and ELISA, suitable for mitochondrial division, apoptosis, and disease mechanism studies. Reliable performance for research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-DRP1-S616 Rabbit pAb
Remarks/Alias	DLP1; DRP1; DVLP; EMPF; OPA5; EMPF1; DYMPLE; HDYNIV; Phospho-DRP1-S616
Species	Human
GeneID (Human)	10059
GeneID	10059
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around S616 of human DRP1 (NP_036192.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O00429
Protein Weight	82 kDa
Shipping	Ice bag