

Anti-Phospho-Egfr-S1070 Rabbit Polyclonal Antibody For Wb, Elisa - P00533

Item Number: CM0000642



Introduction

Phospho-EGFR-S1070 Rabbit pAb detects human EGFR phosphorylated at serine 1070. Suitable for WB and ELISA, this antibody targets a key phosphorylation site in EGFR signaling, ideal for cancer and cell signaling research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-EGFR-S1070 Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-S1070
Species	Human
Gene ID (Human)	1956
Gene ID	1956
Immunogen	Synthetic peptide A phospho specific peptide corresponding to residues surrounding S1070 of human EGFR
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag