


Anti-Phospho-Egfr (T669) Polyclonal Antibody For Wb, Ip, Elisa - P00533

Item Number: CM0000451



Introduction

Anti-Phospho-EGFR (T669) rabbit polyclonal antibody, validated for Western blot, immunoprecipitation, and ELISA. Specifically detects human EGFR phosphorylated at threonine 669, enabling precise study of EGFR signaling cascades in oncology and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-T669 Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-T669
Species	Human
Gene ID (Human)	1956
Gene ID	1956
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T669 of human EGFR (NP_005219.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag