


Anti-Phospho-Egfr-Y1068 Rabbit Monoclonal Antibody For Wb, Elisa - P00533

Item Number: CM0001967



Introduction

Phospho-EGFR-Y1068 rabbit monoclonal antibody specific for EGFR phosphorylated at tyrosine 1068. Validated for WB, ELISA on human samples. Theoretical MW 134 kDa. Ideal for signal transduction and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y1068 Rabbit mAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1068
Species	Human
Gene ID (Human)	1956
Gene ID	1956
Immunogen	Synthetic Peptide A phospho specific peptide corresponding to residues surrounding Y1068 of human EGFR (P00533).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag