


Anti-Phospho-Egfr-Y1172 Rabbit Pab For Wb, Elisa - P00533

Item Number: CM0001853



Introduction

Phospho-specific rabbit polyclonal antibody targeting Y1172 of human EGFR (134 kDa). Validated for WB and ELISA. Ideal for studying EGFR signaling in cancer and cellular responses.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y1172 Rabbit pAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1172
Species	Human
Gene ID (Human)	1956
Gene ID	1956
Immunogen	Synthetic peptide A phospho specific peptide corresponding to residues surrounding Y1172 of human EGFR
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag