


Anti-Phospho-Egfr-Y845 Monoclonal Antibody For Wb, If/Icc, Elisa - P00533

Item Number: CM0001238



Introduction

Rabbit monoclonal antibody specific for EGFR phosphorylated at tyrosine 845. Validated for WB, IF/ICC, and ELISA; ideal for investigating EGFR signaling pathways in oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y845 Rabbit mAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around Y845 of human EGFR (NP_005219.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag