

Anti-Phospho-Eif2A-S51 Monoclonal Antibody For Wb, Elisa - P05198

Item Number: CM0001487



Introduction

Rabbit monoclonal antibody specific to phosphorylated serine 51 of human eIF2 α (EIF2S1). Validated for WB and ELISA applications with cross-reactivity to mouse and rat. Ideal for investigating the integrated stress response and translation regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF2 α -S51 Rabbit mAb
Remarks/Alias	EIF2; EIF-2; EIF2A; EIF-2A; EIF-2alpha; Phospho-eIF2 α -S51
Species	Human
GeneID	1965
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05198
Protein Weight	36 kDa
Shipping	Ice bag