


Anti-Phospho-Eif2A (S51) Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P05198

Item Number: CM0001600



Introduction

Rabbit polyclonal antibody detecting endogenous phosphorylated eIF2 α (Ser51) in human, mouse, rat. Validated for WB, IP, ELISA. Ideal for integrated stress response and translation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF2 α -S51 Rabbit pAb
Remarks/Alias	EIF2; EIF-2; EIF2A; EIF-2A; EIF-2alpha; Phospho-eIF2 α -S51
Species	Human
GeneID (Human)	1965
GeneID	1965
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S51 of human eIF2 α (NP_004085.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05198
Protein Weight	36kDa
Shipping	Ice bag