


Anti-Phospho-Eif4B-S422 Polyclonal Antibody For Wb, Elisa - P23588

Item Number: CM0000632



Introduction

Phospho-specific rabbit polyclonal antibody against human EIF4B phosphorylated at S422. Validated for Western blot and ELISA, recognizes human EIF4B only, with a predicted molecular weight of 69 kDa. Ideal for translation initiation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EIF4B-S422 Rabbit pAb
Remarks/Alias	EIF-4B; PRO1843; Phospho-EIF4B-S422
Species	Human
Gene ID (Human)	1975
Gene ID	1975
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S422 of human EIF4B (NP_001408.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P23588
Protein Weight	69kDa
Shipping	Ice bag