

Anti-Phospho-Eif4E-S209 Rabbit Monoclonal Antibody For Wb, Elisa - P06730

Item Number: CM0001016



Introduction

Phospho-specific rabbit monoclonal antibody targeting eIF4E (S209) for Western blot and ELISA. Cross-reacts with human and rat; detects ~25 kDa eukaryotic translation initiation factor 4E, a key cap-binding protein.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF4E-S209 Rabbit mAb
Remarks/Alias	CBP; EIF4F; AUTS19; EIF4E1; eIF-4E; EIF4EL1; Phospho-eIF4E-S209
Species	Human
GeneID (Human)	1977
GeneID	1977
Immunogen	Synthetic Peptide: A phospho specific peptide corresponding to residues surrounding S209 of human Phospho-eIF4E-S209 (NP_001959.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P06730
Protein Weight	25 kDa
Shipping	Ice bag