


Anti-Phospho-Eif4Ebp1-S65 Rabbit Polyclonal Antibody For Wb, Elisa - Q13541

Item Number: CM0000431



Introduction

Phospho-eIF4EBP1-S65 rabbit polyclonal antibody validated for WB and ELISA. Detects human phosphorylated 4E-BP1 at Ser65, a key mTORC1 target. Ideal for translation regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF4EBP1-S65 Rabbit pAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; Phospho-eIF4EBP1-S65
Species	Human
GeneID (Human)	1978
GeneID	1978
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S65 of human Phospho-eIF4EBP1-S65 (NP_004086.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag