

Anti-Phospho-Eif4Ebp1-T36 Polyclonal Antibody For Wb, Elisa - Q13541

Item Number: CM0001342



Introduction

Phospho-specific rabbit polyclonal antibody against eIF4EBP1 Thr36, validated for WB and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying translation regulation and mTOR signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF4EBP1-T36 Rabbit pAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; Phospho-eIF4EBP1-T36
Species	Human
GeneID (Human)	1978
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T36 of human eIF4EBP1 (NP_004086.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag