


Anti-Phospho-Eif4Ebp1-T70 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13541

Item Number: CM0000495



Introduction

Phospho-specific rabbit polyclonal antibody against eIF4EBP1 (Thr70). Validated in WB, IHC-P, ELISA with human, mouse, and rat reactivity. Ideal for studying translation initiation, mTORC1 signaling, and insulin pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eIF4EBP1-T70 Rabbit pAb
Remarks/Alias	BP-1; 4EBP1; 4E-BP1; PHAS-I; Phospho-eIF4EBP1-T70
Species	Human
Human Gene ID	1978
Immunogen	A synthetic phosphorylated peptide around T70 of human Phospho-eIF4EBP1-T70 (NP_004086.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13541
Protein Weight	13kDa
Shipping	Ice bag