

Anti-Phospho-Enos-S1177 Polyclonal Antibody For Wb, Ihc-P, Elisa - P29474

Item Number: CM0000673



Introduction

Rabbit polyclonal antibody against human eNOS phosphorylated at Ser1177. Suitable for WB, IHC-P, and ELISA, detecting human and mouse eNOS. eNOS produces nitric oxide critical for vascular relaxation and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eNOS-S1177 Rabbit pAb
Remarks/Alias	eNOS; ECNOS; Phospho-eNOS-S1177
Species	Human
Gene ID (Human)	4846
Gene ID	4846
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S1177 of human eNOS (NP_000594.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P29474
Protein Weight	133kDa
Shipping	Ice bag