

Anti-Phospho-Enos (S1177) Rabbit Monoclonal Antibody For Wb, Elisa - P29474

Item Number: CM0001962



Introduction

Specific for the phosphorylated Ser1177 residue of human eNOS, this rabbit monoclonal antibody is validated for Western blot and ELISA, supporting research in endothelial function and cardiovascular disease. SWISS P29474.

[Learn More](#)

Parameter	Value
Product Name	Phospho-eNOS-S1177 Rabbit mAb
Remarks/Alias	eNOS; ECNOS; Phospho-eNOS-S1177
Species	Human
GeneID (Human)	4846
GeneID	4846
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S1177 of human eNOS (NP_000594.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P29474
Protein Weight	133 kDa
Shipping	Ice bag