


## Anti-Phospho-ErbB3/Her3-Y1222 Monoclonal Antibody For Wb, Elisa - P21860

Item Number: CM0001760



### Introduction

Phospho-ErbB3/HER3-Y1222 rabbit monoclonal antibody validated for WB and ELISA, reacting with human ErbB3. Targets phosphorylation at Y1222, ideal for phospho-specific detection in HER3 signaling studies.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-ErbB3/HER3-Y1222 Rabbit mAb                                                                                                   |
| Remarks/Alias    | HER3; FERLK; LCCS2; VSCN1; ErbB-3; c-erbB3; erbB3-S; MDA-BF-1; c-erbB-3; p180-ErbB3; p45-sErbB3; p85-sErbB3; Phospho-ErbB3/HER3-Y1222 |
| Species          | Human                                                                                                                                 |
| GeneID (human)   | 2065                                                                                                                                  |
| GeneID           | 2065                                                                                                                                  |
| Immunogen        | Synthetic Peptide - A synthetic phosphorylated peptide around Y1222 of human ERBB3 (P21860).                                          |
| Source           | Rabbit                                                                                                                                |
| Category         | Phosphorylated Antibodies                                                                                                             |
| Application      | WB, ELISA                                                                                                                             |
| Cross Reactivity | Human                                                                                                                                 |
| SWISS            | P21860                                                                                                                                |
| Protein Weight   | 148kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                               |