

## Anti-Phospho-ErbB3/Her3-Y1289 Rabbit Monoclonal Antibody For Wb, Elisa - P21860

Item Number: CM0001055



### Introduction

Phospho-specific rabbit monoclonal antibody against ErbB3 (HER3) Y1289, validated for WB and ELISA. Detects human ErbB3 phosphorylated at tyrosine 1289, suitable for studying receptor tyrosine kinase signaling in epithelial and cancer research.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-ErbB3/HER3-Y1289 Rabbit mAb                                                                                                   |
| Remarks/Alias    | HER3; FERLK; LCCS2; VSCN1; ErbB-3; c-erbB3; erbB3-S; MDA-BF-1; c-erbB-3; p180-ErbB3; p45-sErbB3; p85-sErbB3; Phospho-ErbB3/HER3-Y1289 |
| Species          | Human                                                                                                                                 |
| GenelD           | 2065                                                                                                                                  |
| Immunogen        | A synthetic phosphorylated peptide around Y1289 of human ErbB3                                                                        |
| Source           | Rabbit                                                                                                                                |
| Category         | Phosphorylated Antibodies                                                                                                             |
| Application      | WB, ELISA                                                                                                                             |
| Cross Reactivity | Human                                                                                                                                 |
| SWISS            | P21860                                                                                                                                |
| Protein Weight   | 148kDa                                                                                                                                |
| Shipping         | Ice bag                                                                                                                               |